

Micro Urine Collection

for Babies and Toddlers Not Toilet-Trained

Equipment

- 1x 70ml Yellow-top Container
- 1x Paediatric Urine Collector
- Sterile Water Wipe

Procedure

For children not toilet-trained

- The parent/carer must wash hands thoroughly.
- Cleanse child's genital area thoroughly with sterile water wipe.
- Put on the paediatric urine collector (remove backing paper from adhesive by lifting the lower left or right hand corner).
- Remove centre piece to expose the opening.
- For males, attach the bag allowing the penis to project through the opening.
- For females, attach the bag so the opening covers the upper part of the genitalia.
- As soon as possible after urine has been passed (minimum requirement is 1ml), transfer to the sterile container.
- Remove lid and hold bag over the container.
- Empty the bag by cutting the projection at the bottom corner and allowing the urine to flow into the container. Immediately replace the lid, ensuring correct threading and secure closure.
- **Write the patient's surname, first name, date of birth, and the word 'Urine' on the container. Also write the time and date the sample was collected.**
- Store in the refrigerator in a sealed bag and deliver to the collection centre within 24 hours of collection.